

I-TECH
MEDICAL DIVISION

**Our mission, your
freedom:
- pain
+ rehabilitation**



Distributed by:

InnoMed Spring Ltd.

Tel: +852 3490 4862

Email: info@innomedspring.com.hk

© 2026 InnoMed Spring Ltd. | All rights reserved.






IACER, the company.

IACER is an Italian company with over 50 years of history. Since 2004, through the I-Tech Medical Division brand, we have been **committed to research, innovation and home-patient care.**

We design and manufacture **certified medical devices for home rehabilitation**—including PEMF therapy, electrotherapy, ultrasound therapy, pressotherapy, and diagnostics, as well as advanced equipment for physiotherapists, such as diathermy.

At I-Tech Medical Division we **support people with musculoskeletal and vascular conditions throughout their entire therapeutic journey.** Our goal is to help them recover fully, with continuous support also provided via our **telerehabilitation app.**

Choose the solution that suits you best and discover our medical devices



MAGNETOTHERAPY / PEMF

Reduce pain and inflammation with magnetotherapy. The ideal ally for treating traumatic and chronic conditions affecting bones, muscles and joints. The action of pulsed electromagnetic fields will help your body's cells restore their balance and proper functioning.

RECOMMENDED DURATION OF TREATMENT



At least 2-4 hours per day, depending on the pathology.



45-90 days.
Repeat, if necessary.

THERAPEUTIC INDICATIONS

INFLAMMATORY PATHOLOGIES - OSTEOPOROSIS - ARTHROSIS - ARTHRITIS
INFLAMMATION - BONE OEDEMA - FRACTURES - LUMBAGO - ALGODYSTROPHY
JOINT PAIN - ULCERS - NEUROPATHIES

LAMAGNETO PRO

Even simpler, even more secure. Lock the program with a password and always start the correct therapy.

63 PRESET PROGRAMS

to make your therapy easy

CUSTOMISED TREATMENTS

thanks to 10 free memories

2 TREATMENTS PER SESSION

and get more results in less time

LOCK WITH PASSWORD

and you are sure you are not varying the correct therapy parameters



TRACK THE THERAPIES PERFORMED

with integrated counter

INTENSITY UP TO 200 GAUSS

for each of the 2 channels

FREQUENCY FROM 1 TO 120 HZ

to set precisely the parameters prescribed by your doctor



MD Medical Device



Connect it to the APP



Recommended by doctors.

Included



PROGRAMS (LAMAGNETO PRO AND LAMAGNETO X)

01 Osteoporosis	09 Epicondylitis	17 Coxarthrosis	25 Tendinitis
02 Arthrosis	10 Epitrocleitis	18 Muscle atrophies	26 Chondropathy
03 Arthritis	11 Pseudoarthrosis	19 Muscle treatment	27 Anti-inflammatory
04 Cervical arthrosis	12 Lumbago	20 Osteonecrosis	28 Whiplash
05 Joint pain	13 Lumbosacral pain	21 Algodystrophy	29 Scarring
06 Cervicalgia	14 Shoulder arthrosis	22 Cartilage lesion	30 Skin ulcers
07 Chronic pain	15 Knee arthrosis	23 Ligament injury	31 Discopathy
08 Fractures	16 Scapulohumeral periarthritis	24 Bone oedema	32 Myalgia



CV2S2000 SOLENOID PAIR

Pair of professional solenoids for high-intensity therapy (solenoid diameter: 11 cm)

Compatible with LaMagneto, LaMagneto Pro, LaMagneto X, Mag2000, Mag2000 Premium, Mag2000 Plus



CE MD

F3S2000 THERAPEUTIC BELT

Therapeutic elastic belt for low-frequency magnetotherapy with 3-solenoid cable.

Compatible with LaMagneto, LaMagneto Pro, LaMagneto X, Mag2000, Mag2000 Premium, Mag2000 Plus.



CE MD

THERAPEUTIC BELT F2S700

Elastic therapeutic belt for low-frequency magnetotherapy with 2-solenoid cable.

Compatible with Mag700.



CE MD

TAP2000

The 40x40cm mat with 6 solenoids, ideal for your back treatments.

Compatible with LaMagneto, LaMagneto Pro, LaMagneto X, Mag2000, Mag2000 Premium, Mag2000 Plus.



CE MD

CIL-MAG

The cylindrical professional solenoid which allows you to comfortably treat the lower limbs.

Compatible con Mag2000, Mag2000 Premium, Mag2000 Plus.



CE MD

ORTHOMAG

Are you always on the move?

With Orthomag you have therapy wherever you want, in total freedom.

**PORTABLE,
WITH RECHARGEABLE
BATTERY**

for therapy wherever you want, whenever you want



**INTENSITY
UP TO 25 GAUSS**

**FLEXIBLE
APPLICATOR**



CE 0197

Included



TREATMENTS

PROGRAM L (50 HZ)

Algodystrophies
Gonarthrosis
Cartilage degeneration
Fractures
Cartilage injuries
Osteoarthrosis
Osteonecrosis
Osteoporosis
Delayed calcification
Pseudo=arthrosis

PROGRAM H (75 HZ)

Scar adherence
Arthritis
Arthrosis
Bursitis
Brachialgy
Capsulitis
Cervical pain
Whiplash
Chondropathy
Contusions
Coxarthrosis
Articular pain
Rheumatic disorders
Back pain
Epicondylitis
Epitrochleitis
Discal Hernia
Plantar fasciitis
Lumbago
Meniscopathy
Metatarsal pain
Periarthritis
Pubalgy
Rheumatisms
Rhizarthrosis
Sciatalgy
Inflammatory diseases
Muscle strains

ELECTROTHERAPY

Recover faster, reducing pain.

Electrotherapy transmits micro-impulses to the human body, varying in frequency pulse width and intensity, depending on the treatment to carry out.

RECOMMENDED TREATMENT DURATION



To find out the duration of therapy we recommend that you consult your doctor and follow the instructions in our user manual.

THERAPEUTIC INDICATIONS

MUSCLE RECOVERY – PAIN REDUCTION – INCONTINENCE TREATMENT – TREATMENT OF DENERVATED MUSCLE – FLUID DRAINAGE

MIO-CARE TENS

Reduce pain quickly.
In case of both chronic and acute pain.

PORTABLE
BATTERY-
POWERED



CE 0197

20 TENS programs

-

-

2 independent channels


Compensated biphasic square wave

Belt clip

Maximum intensity: 100mA

Internal rechargeable Ni-Mh battery

Included



MIO-CARE FITNESS

Find pain relief, reduce blemishes and strengthen your muscles. All with Mio-Care Fitness.

PORTABLE
BATTERY-
POWERED



CE 0197

20 TENS programs

20 FITNESS programs

15 BEAUTY programmes

2 independent channels


Compensated biphasic square wave

Belt clip

Maximum intensity: 100mA

Internal rechargeable Ni-Mh battery

Included



I-TECH PHYSIO

The complete and programmable device to treat pain and rehabilitate your muscles.



Recommended by doctors.



Clinical studies on prehabilitation prior to arthroplasty

PROGRAMS

TENS

- 01 Rapid tens
- 02 Endorphinic tens
- 03 Tens at maximum values
- 04 Anti-inflammatory
- 05 Cervicalgia/Myotensive headache
- 06 Lumbago/Sciatica
- 07 Sprains/Contusions
- 08 Vascularisation
- 09 Decontracting
- 10 Hand and wrist pain
- 11 Plantar stimulation
- 12 Epicondylitis
- 13 Epitrocleitis
- 14 Periarthritis

REHA

- 01 Iontophoresis L
- 02 Iontophoresis M
- 03 Iontophoresis H
- 04 Microcurrent
- 05 Haematoma
- 06 Edema
- 07 Modulated Tens
- 08 Burst Tens
- 09 Atrophy prevention
- 10 Atrophy
- 11 Denervated quadra 50 ms
- 12 Denervated quadra 100 ms
- 13 Denervated quadra 150 ms
- 14 Denervated quadra 200 ms
- 15 Denervated quadra 250 ms
- 16 Denervated triangular 50 ms
- 17 Denervated triangular 100 ms
- 18 Denervated triangular 150 ms
- 19 Denervated triangular 200 ms
- 20 Denervated triangular 250 ms
- 21 Interferential
- 22 Tens with amplitude modulation
- 23 Alternating Tens
- 24 Tens with frequency modulation 1
- 25 Tens with frequency modulation 2
- 26 Tens with amplitude modulation
- 27 Kotz



2 INDEPENDENT CHANNELS

easily programmable

DENERVATED MUSCLE, IONTOPHORESIS INCONTINENCE

INTERFERENTIAL AND KOTZ WAVE

in addition to the classic compensated biphasic square, monophasic square and triangular



Connect it to the APP

BELT CLIP
for greater comfort

MAXIMUM INTENSITY 100 mA

PORTABLE BATTERY-POWERED



Included

Possibility to add

Anal probe
Vaginal probe



NEMS

- 01 Warm-up
- 02 Upper limb resistance and trunk
- 03 Lower limb Endurance
- 04 Resistance upper limbs and trunk
- 05 Lower limb strength
- 06 Basic strength upper limbs and trunk
- 07 Basic strength lower limbs
- 08 Fast strength upper limbs and trunk
- 09 Fast strength lower limbs
- 10 Explore strength upper limbs and trunk
- 11 Explosive strength lower limbs
- 12 Deep capillarisation
- 13 Muscle recovery
- 14 Agonist-antagonist
- 15 Sequential tonic contractions Upper limbs and trunk

- 16 Tonic contractions sequential lower limbs contractions
- 17 Phasic contractions sequential upper limbs and trunk
- 18 Phasic contractions sequential lower limbs and trunk
- 19 Decontracting
- 20 Deep tissue massage
- 21 EMS

BEAUTY

- 01 Firming upper limbs and trunk
- 02 Lower limb toning
- 03 Upper limb toning and trunk
- 04 Lower limb toning
- 05 Upper limb definition and trunk
- 06 Lower limb definition
- 07 Modelling
- 08 Microlifting
- 09 Abdomen lipolysis
- 10 Lipolysis thighs
- 11 Lipolysis buttocks and hips
- 12 Lipolysis arms
- 13 Tissue elasticity
- 14 Capillarisation

INCONTINENCE

- 01 Stress incontinence*
- 02 Stress incontinence*
- 03 Stress incontinence*
- 04 Urge incontinence*
- 05 Urge incontinence*
- 06 Urge incontinence*
- 07 Mixed incontinence*
- 08 Mixed incontinence*
- 09 Mixed incontinence*

*Probes not included

+ 13 FREE MEMORIES

ULTRASOUND THERAPY

With ultrasound you can find relief from pain in just a few sessions.

By transferring energy to the tissues of our body, they promote microcirculation and help eliminate excess subcutaneous fluid properly.

RECOMMENDED TREATMENT DURATION



5 to 15 minutes
per day.



10 days.
Potentially
repeatable.

THERAPEUTIC INDICATIONS

MYOFASCIAL PAIN – MUSCLE PAIN – CARPAL TUNNEL – VENOUS ULCERS –
PRESSURE ULCERS – LUMBAR PAIN – ARTHROSIS – EPICONDYLITIS – EPITROCH-
LEITIS – TENDINITIS – CUBITAL TUNNEL – LUMBAR STENOSIS – SCIATICA

POWERSONIC

Do therapy where you want, when you want. And reduce pain and inflammation in just a few sessions.



MAXIMUM INTENSITY
1.6 W/CM²
plugged into the socket

PORTABLE
BATTERY-POWERED

1 MHZ ULTRASOUND
like a professional treatment

ADJUSTABLE
DURATION
AND INTENSITY
for 9 different combinations



Connect it to the **APP**



Recommended
by doctors.

CE 0068

TREATMENTS

Myofascial pain	Pressure ulcers	Epitrocleitis	Sciatica
Muscle pain	Lumbar pain	Tendinitis	
Carpal tunnel	Arthrosis	Cubital tunnel	
Venous ulcers	Epicondylitis	Lumbar stenosis	

I-TECH UT²

Doublehandpiece(5cm² and 1cm²) and choice of working frequency (1-3MHz) to adapt the treatment to the patient's therapeutic needs.



Included

Add Troll-M!
3-shelf trolley with 4 swivel castors and 1 socket power socket.



CE 0068

I-TECH UE

A complete device. Ultrasound therapy, electrotherapy with the possibility of combining the two therapies for even greater benefits.



Included

Add Troll-M!
3-shelf trolley with 4 swivel castors and 1 socket power socket.



CE 0068

PRESSOTHERAPY

Swollen or stiff legs?

Pressotherapy will help you reduce swelling, limb stiffness and remove toxins.

It also promotes proper venous and lymphatic return in arms and legs.

RECOMMENDED TREATMENT DURATION



From 20 to 60
minutes per day



30 days

THERAPEUTIC INDICATIONS

EDEMA – LYMPHEDEMA – VENOUS ULCERS – VENOUS INSUFFICIENCY

MUSCLE RECOVERY

HOW DOES IT WORK?

Pressotherapy exerts sequential compression on the limbs and abdomen with a disto-proximal movement, i.e. from the periphery towards the centre of the body. This action favours the movement of blood and lymph in the presence of stagnation or circulation difficulties.

I-PRESS

I-Press helps you improve venous and lymphatic circulation of the limbs and treat certain pathologies of the circulatory system.

4 COMPRESSING CHAMBERS

overlapped which inflate in sequence from foot to thigh

CHOOSE CONFIGURATION AND ACCESSORIES

customise the treatment as you wish

NEW QUIETER COMPRESSOR

EASY AND INTUITIVE

set session duration and pressure and you are ready for the therapy



Clinical studies on edema and recovery post-arthroplasty with our device



Connect it to the **APP**



Recommended by doctors.



MD Medical Device

CONFIGURATIONS

I-PRESS LEG1

- Machine body
- n. 1 legging
- n. 1 4-chamber 4-ending tube
- n. 1 connector
- n. 1 pressure orthotic



I-PRESS LEG2

- Machine body
- n. 2 leggings
- n. 2 4-chamber 4-ending tubes
- n. 1 splitter
- n. 1 connector
- n. 1 pressure orthotic



I-PRESS ARMI

- Machine body
- n. 1 arm cuff
- n. 1 4-chamber 4-ending tube
- n. 1 connector



All 3 configurations are available in sizes M and L and can be further customised with the full range of pressotherapy accessories.

TREATMENTS

Pressure (mmHg)

Time (minutes)

	Pressure (mmHg)	Time (minutes)
Edema	50	50
Lymphedema	50	30
Venous ulcers	50	60
Venous insufficiency	50	60
Muscle recovery	85	20

LEGGINGS

LEGGING

THIGH CIRCUMFERENCE

ANKLE CIRCUMFERENCE

LEG LENGTH

(measurement taken from inner thigh to heel)

SIZE **M**



Up to 64 cm

Up to 50 cm

Up to 82 cm

SIZE **L**



Up to 70 cm

Up to 50 cm

Up to 98 cm

LEGGING EXTENDERS

Pair of leg extenders



TNT PROTECTION

Pair of hygienic TNT leg protections



PRESSURE ORTHOTIC



ARMBAND

4-compression chamber arm cuff

ARMBAND

ARMPIT CIRCUMFERENCE

WRIST CIRCUMFERENCE

SIZE **M**



Up to 47 cm

Up to 32 cm

SIZE **L**




Up to 56 cm

Up to 40 cm

REATHERM

Diathermy with the neuromodulation effect:
effective, safe, precise.

Optimum results from the very first session.
Whether for rehabilitation, sport, beauty.



PROGRAMS

REHAB

01 Phantom limb	11 TMJ Dysfunction	21 ACL
02 Knee arthrosis	12 Sprain	22 Lymphedema
03 Arthrosynovitis	13 Muscle Distraction	23 Low back pain
04 Bursitis	14 Joint pain	24 Sciatica
05 Neck pain	15 Edema	25 Meniscopathy (meniscal synovitis)
06 Whiplash syndrome	16 Enthesitis	26 Diabetic foot
07 Chondropathy	17 Epicondylitis	27 Piriformis syndrome
08 Muscle contracture	18 Epitrochleitis	28 Muscle tear
09 Contusion	19 Disc hernia	29 Acute tendinitis
10 Muscular sprain	20 Scarring	30 Chronic tendinitis

DIATHERMY

MODULATION 1 ÷ 100 HZ

Comfortable treatment, also in athermia, neuromodulation effect. Not to be confused with the classic duty-cycle.

401 FREQUENCIES (300 ÷ 700 KHZ)

to adapt the treatment to the needs of the individual patient

10 FREE MEMORIES MODE CUSTOM

customised treatment, up to 60 minutes

UNIPOLAR AND BIPOLAR HANDPIECES

ergonomic

50 PRESET PROGRAMS

for areal guide to treatments

PROFESSIONAL TRAINING INCLUDED

LIGHTWEIGHT, POWERFUL, TRANSPORTABLE

3 kg, 200 W

Included



Add Troll-M!

3-shelf trolley with 4 swivel castors and 1 socket power socket.



CE 0068

SPORT

- 31 Rotator cuff tendinopathy
- 32 Ulcer
- 33 Effusion
- 34 Arthroplasty (recovery)
- 35 Fractures

- 01 Muscle contracture
- 02 Contusion
- 03 Muscle recovery
- 04 Sprain
- 05 Hematoma
- 06 Tendon injury
- 07 Superficial massage
- 08 Deep massage
- 09 Muscle sprain
- 10 Muscle tear

BEAUTY

- 01 Cellulite
- 02 Collagen
- 03 Edema
- 04 Draining massage
- 05 Microlifting

+ 10 FREE MEMORIES



I-TECH
Life

Get the most
from your therapy.

**Download
our new app
and start your own
therapeutic journey.**



Manufactured by:

I-TECH
MEDICAL DIVISION

a **IACER** brand

I.A.C.E.R. Srl

via Enzo Ferrari 2

30037 Scorzè (VE) - Italy

Tel.: (+39) 041/5401356

itechmedicaldivision.com

C.F. e P.Iva: IT00185480274

customercare@iacer.it



Distributed by:

InnoMed Spring Ltd.

Tel: +852 3490 4862

Email: info@innomedspring.com.hk

© 2026 InnoMed Spring Ltd. | All rights reserved.